


CMV promoter 1–602
multiple cloning site 651–758
SV40 polyA 811–1194
f1 origin 1332–1638
bla promoter 1692–1816
SV40 promoter 1836–2174
neomycin/kanamycin resistance ORF 2209–3000
HSV-TK polyA 3001–3450
pUC origin 3588–4255



pCMV-Script EX Multiple Cloning Site Region
(sequence shown 620–799)

T3 promoter
 AATTAACCCCTCACTAAAGGGAACAAAAGCTGGAGCTCCACCGCGGTGGCGGCCGCTCTA...

Sac I BstX I Sac II Not I

Srf I BamH I Pst I EcoR I EcoR V Hind III Acc I/Sal I

...GCCCGGGCGGATCCCCGGGCTGCAGGAATTCGATATCAAGCTTATCGATACCGTCGAC...

Xho I Apa I Kpn I

...CTCGAGGGGGGGCCCGGTACCAGGTAAGTGTACCCAATTCGCCCTATAGTGAGTCGTATTAC

T7 promoter

