


CMV promoter 1–602
multiple cloning site 651–758
SV40 polyA 811–1194
f1 origin 1332–1638
bla promoter 1663–1787
SV40 promoter 1807–2145
neomycin/kanamycin resistance ORF 2180–2971
HSV-TK polyA 2972–3421
pUC origin 3559–4226



pCMV-Script Multiple Cloning Site Region
 (sequence shown 620–799)

T3 promoter
AATTAACCCCTCACTAAAGGGAACAAAAGCTGGAGCTCCACCGCGGTGGCGGCCGCTCTA...

Sac I BstX I Sac II Not I

Srf I BamH I Pst I EcoR I EcoR V Hind III Acc I/Sal I

...GCCCGGGCGGATCCCCGGGCTGCAGGAATTCGATATCAAGCTTATCGATACCGTCGAC...

Xho I Apa I Kpn I

...CTCGAGGGGGGGCCCGGTACCAGGTAAGTGTACCCAATTCGCCCTATAGTGAGTCGTATTAC

T7 promoter