


CMV promoter 1–602
SV40 pA 811–1194
f1 origin 1363–1669
bla promoter 1694–1818
SV40 promoter 1838–2176
neomycin/kanamycin resistance ORF 2211–3005
HSV-TK polyA 2962–3420
pUC origin 3590–4257

Note *The pCMV-PCR vector provided with Catalog #200480 has been linearized and LIC-adapted, (LIC overhangs are indicated). The sequence shown in gray is absent from the LIC-adapted vector. The positions listed above correspond to the pCMV-PCR vector **prior** to LIC adaptation.*



pCMV-PCR Multiple Cloning Site Region
(sequence shown 620–830)

