



**Comments for pCMV/myc/cyto/GFP**  
5611 nucleotides

- CMV immediate-early promoter: bases 12-627
- CMV priming site: bases 544-564
- GFP ORF: bases 667-1393
- myc* epitope: bases 1403-1432
- BGH Reverse priming site: bases 1468-1485
- BGH polyadenylation sequence: bases 1467-1681
- f1 origin: bases 1744-2157
- SV40 promoter/origin: bases 2222-2533
- Neomycin (G418) resistance gene (ORF): bases 2541-3335
- SV40 polyadenylation sequence: bases 3351-3590
- pUC origin: bases 3775-4448
- Ampicillin resistance gene (ORF): 4593-5453